

Cylinder



Hydrometer Jar, Round base, Plain

- Without calibration markup
- Without Spout

Item Code	Capacity (ml)	Pack of
138.271.02	250	2
138.271.03	500	2
138.271.04	1000	1

Hydrometer Jar, Round base

- With circle mark on the top
- Printed with capacity & 20°C
- Without spout
- With rim

Item Code	Capacity (ml)	Pack of
138.270.01	1000	1

3.3
BORO
GLASS

BLUE
PRINTING

Chemical
Resistant

A
CLASS



ENQUIRE



ENQUIRE



Nessler Cylinder, for colour comparison Class-A

Item Code	Capacity (ml)	Sub Division (ml)	Tolerance (± ml)	Pack of
143.202.01	50	25 & 50	0.4	10
143.202.02	100	50 & 100	0.8	10

3.3
BORO
GLASS

BLUE
PRINTING

Chemical
Resistant

A
CLASS



ENQUIRE