


# Volumetric Flask



ENQUIRE

## Volumetric Flask, PE Stopper, Class A

- Complies with USP Standard
- With Individual work Certificate and Serial Number
- QR Coded Flask
- With PE Stopper

Item Code	Capacity (ml)	N/S	Tolerance (±)	Pack of
QR.130.220.01	5	10/19	0.020	1
QR.130.220.02	10	10/19	0.020	1
QR.130.220.02A	20	10/19	0.030	1
QR.130.220.03	25	10/19	0.030	1
QR.130.220.04	50	12/21	0.050	1
QR.130.220.04B	50	14/23	0.050	1
QR.130.220.05B	100	12/21	0.080	1
QR.130.220.05	100	14/23	0.080	1
QR.130.220.06	200	14/23	0.100	1
QR.130.220.07	250	14/23	0.120	1
QR.130.220.08	500	19/26	0.150	1
QR.130.220.09	1000	24/29	0.300	1
QR.130.220.10	2000	29/32	0.500	1

3.3  
BORO  
GLASS

Chemical  
Resistant

AMBER  
PRINTING

A  
CLASS

PE

AMBER  
RING

INDIVIDUAL  
CERTIFICATE



ENQUIRE

## Volumetric Flask, PE Stopper, Class A

- Complies with USP Standard
- With Batch Certificate
- With PE Stopper

Item Code	Capacity (ml)	N/S	Tolerance (±)	Pack of
130.461.01	5	10/19	0.020	2
130.461.02	10	10/19	0.020	2
130.461.03	20	10/19	0.030	2
130.461.04	25	10/19	0.030	2
130.461.05	50	12/21	0.050	2
130.461.06	50	14/23	0.050	2
130.461.07	100	12/21	0.080	2
130.461.08	100	14/23	0.080	2
130.461.09	200	14/23	0.100	2
130.461.10	250	14/23	0.120	2
130.461.11	500	19/26	0.150	2
130.461.12	1000	24/29	0.300	2
130.461.13	2000	29/32	0.500	1

3.3  
BORO  
GLASS

Chemical  
Resistant

AMBER  
PRINTING

A  
CLASS

PE

AMBER  
RING

BATCH  
CERTIFICATE