



All Glass Filter Holder

For 90 mm Disc Filters

Applications :

- Used in vacuum filter aqueous, organic or corrosive liquids for particulate contamination analysis
- Recommended for HPLC solvent filtration.

Specifications:

Materials	Borosilicate glass funnel, base and tubulated cap; PTFE-coated stainless steel screen; anodized aluminum spring clamp
Filter Diameter, mm	90
Funnel Capacity, L	1
Outlet Fitting	6 mm (1/4 in.) O.D. tubulated cap sidearm to vacuum

Ordering Information:

Cat. No.	Description
263.202.01	Glass Filter Holder with PTFE Coated Stainless steel screen, 90 mm

Replacement Parts

Cat. No.	Description
263.245.01	Funnel, 1L, 90 mm, ground glass seal
263.245.02	Vacuum Base and Cap, 90 mm
263.245.03	Stainless Steel Screen, PTFE coated, 90 mm
263.245.04	Gasket, PTFE 90 mm
263.245.05	Spring Clamp, 90 mm anodized aluminum



Cat. No.	Description
263.202.02	Glass Filter Holder with frit ted disc

Replacement Parts

Cat. No.	Description
263.245.01	Funnel, 1L, 90 mm, ground glass seal
263.248.02	Cap glass base with fritted disc
263.245.05	Spring Clamp, 90 mm anodized aluminum

Optional Accessories

Cat. No.	Description
258.245.04	Ground Joint Flask, 1 L
263.245.06	Bottle Vacuum 2 L
263.245.07	Bottle Vacuum 5 L
258.245.08	Ground joint Flask, 4 L
258.245.09	Ground Joint Flask 5 L



ENQUIRE