

Adapter

- Our wide range of interchangeable jointed glassware adapters ensure a secure and airtight connection for efficiency in your apparatus.
- All adapters, ranging between reduction, expansion, multi neck, inlet, short & long stem with/without drip tips in various sizes and shapes, are manufactured from ASTM E438, TYPE 1: Class A Boro 3.3 glass that meets DIN/ISO 3585 & USP Type I standards.



Adapter, Bushing with Drip Tip

- Thick wall and reinforced rim.
- Full length (larger standard taper size) outside of the bushing and shorter length inside of the bushing
- To compactly and conveniently assemble dissimilar sized joints glassware

| Item Code | Socket Size | Cone Size | Pack of |
|------------|-------------|-----------|---------|
| 013.230.00 | 10/19 | 14/23 | 10 |
| 013.230.01 | 14/23 | 19/26 | 10 |
| 013.230.02 | 19/26 | 29/32 | 10 |

ASTM Sizes

| Item Code | Socket Size | Cone Size | Pack of |
|-------------|-------------|-----------|---------|
| 013.478.01 | 14/20 | 19/22 | 1 |
| 013.478.01A | 14/20 | 24/40 | 1 |
| 013.478.02 | 14/35 | 24/40 | 1 |


3.3
BORO
GLASS



ENQUIRE



ENQUIRE



Adapter, Reduction

- Complies as per DIN 12257
- Small sized socket at the top and larger sized cone at the base to connect dissimilar size ground glass joints

| Item Code | Socket Size | Cone Size | Pack of |
|-------------|-------------|-----------|---------|
| 013.202.01* | 14/23 | 14/23 | 10 |
| 013.202.01A | 14/23 | 19/26 | 10 |
| 013.202.02* | 14/23 | 24/29 | 10 |
| 013.202.03 | 14/23 | 29/32 | 10 |
| 013.202.04* | 19/26 | 24/29 | 10 |
| 013.202.05 | 19/26 | 29/32 | 10 |
| 013.202.06* | 19/26 | 34/35 | 10 |
| 013.202.07* | 19/26 | 40/38 | 10 |
| 013.202.08* | 19/26 | 45/40 | 10 |
| 013.202.09* | 19/26 | 55/44 | 10 |
| 013.202.10* | 24/29 | 29/32 | 10 |
| 013.202.11* | 24/29 | 34/35 | 10 |
| 013.202.12* | 24/29 | 40/38 | 10 |
| 013.202.13* | 24/29 | 45/40 | 10 |
| 013.202.14* | 24/29 | 50/42 | 10 |
| 013.202.15* | 24/29 | 55/44 | 10 |
| 013.202.16* | 29/32 | 34/35 | 10 |
| 013.202.17* | 29/32 | 40/38 | 10 |
| 013.202.18 | 29/32 | 45/40 | 10 |
| 013.202.19* | 34/35 | 40/38 | 10 |
| 013.202.20* | 34/35 | 45/40 | 10 |
| 013.202.21* | 34/35 | 50/42 | 10 |

ASTM Sizes

| Item Code | Socket Size | Cone Size | Pack of |
|-------------|-------------|-----------|---------|
| 013.470.02A | 14/20 | 19/22 | 1 |
| 013.470.02 | 14/20 | 24/40 | 1 |
| 013.470.04 | 19/22 | 24/40 | 1 |
| 013.470.10 | 24/40 | 29/42 | 1 |
| 013.470.13 | 24/40 | 45/50 | 1 |

3.3
BORO
GLASS



1 x
SOCKET

1 x
CONE



ENQUIRE



ENQUIRE

*Non DIN/ISO