

Funnel, Separating, Pear (ISO) Shape With Glass Stopcock

- Complies with ISO 4800 Standard
- With PE Stopper


Item Code	Capacity (ml)	N/S	Stopcock Type	Bore Size	Pack of
150.202.02	50	19/26	Solid Glass	2.5	1
150.202.03	100	19/26	Solid Glass	2.5	1
150.204.03A*	125	19/26	Solid Glass	2.5	1
150.202.04	250	29/32	Solid Glass	4	1
150.202.05	500	29/32	Solid Glass	4	1
150.202.06	1000	29/32	Solid Glass	6	1
150.202.07	2000	29/32	Solid Glass	6	1

3.3
BORO
GLASS



ENQUIRE

* Non DIN



Funnel, Separating, Pear (ISO) Shape With Needle Valve Stopcock

- Complies with ISO 4800 Standard
- With PE Stopper
- Bore 0-3 mm

Item Code	Capacity (ml)	N/S	Stopcock Type	Pack of
151.202.02	50	19/26	Needle Valve	1
151.202.03	100	19/26	Needle Valve	1
151.204.03A*	125	19/26	Needle Valve	1
151.202.04	250	29/32	Needle Valve	1
151.202.05	500	29/32	Needle Valve	1
151.202.06	1000	29/32	Needle Valve	1
151.202.07	2000	29/32	Needle Valve	1

3.3
BORO
GLASS



ENQUIRE

* Non DIN



Funnel, Separating, Pear (ISO) Shape With PTFE Stopcock

- Complies with ISO 4800 Standard
- With PE Stopper

Item Code	Capacity (ml)	N/S	Bore (mm)	Stopcock Type	Pack of
152.202.02	50	19/26	2.5	PTFE	1
152.202.03	100	19/26	2.5	PTFE	1
152.202.03A*	125	19/26	2.5	PTFE	1
152.202.04	250	29/32	4	PTFE	1
152.202.05	500	29/32	4	PTFE	1
152.202.06	1000	29/32	6	PTFE	1
152.202.07	2000	29/32	6	PTFE	1

3.3
BORO
GLASS



ENQUIRE

* Non DIN

