

Adapter, Cone for Rubber Tubing

- Right angled connection for the attachment of rubber hose
- With Drip Tip
- With Screw Thread

Item Code	Cone Size	Pack of
036.202.01	14/23	1
036.202.02	19/26	1
036.202.03	24/29	1
036.202.04	29/32	1

ASTM Sizes

Item Code	Cone Size	Pack of
036.475.01S	14/20	1
036.475.03S	24/40	1

3.3
BORO
GLASS



1 x
CONE



ENQUIRE



ENQUIRE



Adapter, Cone with long Stem

- With screw thread, right angle tubing
- Detachable plastic screw thread that makes it compatible for a rubber hose connection

Item Code	Cone Size	Pack of
037.202.01	14/23	1
037.202.02	19/26	1
037.202.03	24/29	1
037.202.04	29/32	1

3.3
BORO
GLASS



1 x
CONE



ENQUIRE



Adapter, Tubing Straight

- With Hooks
- Detachable plastic screw thread that makes it compatible for a rubber hose connection

Item Code	Socket Size	Pack of
038.202.01	14/23	1
038.202.02	19/26	1
038.202.03	24/29	1
038.202.04	29/32	1

ASTM Sizes

- Without Hooks

Item Code	Socket Size	Pack of
020.470.01	14/20	1
020.470.03	24/40	1

3.3
BORO
GLASS



1 x
SOCKET



ENQUIRE



ENQUIRE

