

## Volumetric Flasks, Amber, Class A

- Complies with DIN ISO 1042 Standards
- With Batch certificate
- With PE Stopper

| Item Code    | Capacity (ml) | N/S   | Tolerance (±) | Pack of |
|--------------|---------------|-------|---------------|---------|
| 130.202.01A  | 5             | 10/19 | 0.040         | 2       |
| 130.202.02AB | 10            | 10/19 | 0.040         | 2       |
| 130.202.03A  | 20            | 10/19 | 0.040         | 2       |
| 130.202.04A  | 25            | 10/19 | 0.040         | 2       |
| 130.202.05A  | 50            | 12/21 | 0.060         | 2       |
| 130.202.05AB | 50            | 14/23 | 0.100         | 2       |
| 130.202.06AB | 100           | 12/21 | 0.100         | 2       |
| 130.202.06A  | 100           | 14/23 | 0.100         | 2       |
| 130.202.07A  | 200           | 14/23 | 0.150         | 2       |
| 130.202.08A  | 250           | 14/23 | 0.150         | 2       |
| 130.202.09A  | 500           | 19/26 | 0.250         | 2       |
| 130.202.10A  | 1000          | 24/29 | 0.400         | 2       |
| 130.202.11A  | 2000          | 29/32 | 0.600         | 1       |



ENQUIRE



3.3  
BORO  
GLASS

Chemical  
Resistant

WHITE  
PRINTING

A  
CLASS

PE



BATCH  
CERTIFICATE

## Trapezoidal Volumetric Flask, CLASS A, Amber Glass *New*

- Boro 3.3 glass
- With Batch certificate
- With PE Stopper

| Item Code    | Capacity (ml) | N/S   | Pack of |
|--------------|---------------|-------|---------|
| 135.202.00A  | 1             | 7/16  | 2       |
| 135.202.00BA | 2             | 7/16  | 2       |
| 135.202.01A  | 5             | 10/19 | 2       |
| 135.202.02A  | 10            | 10/19 | 2       |
| 135.202.03A  | 20            | 10/19 | 2       |
| 135.202.04A  | 25            | 10/19 | 2       |
| 135.202.05A  | 50            | 12/21 | 2       |



ENQUIRE



3.3  
BORO  
GLASS

Chemical  
Resistant

WHITE  
PRINTING

A  
CLASS

PE



BATCH  
CERTIFICATE