

# Cylinder



## Mixing Cylinder with Hexagonal base, Class-A

- Complies with DIN EN ISO 4788 standards
- With conformity batch certificate
- With PE stopper

Item Code	Capacity (ml)	N/S	Sub Division (ml)	Tolerance ( $\pm$ ml)	Pack of
142.202.00A*	5	10/19	0.1	0.05	2
142.202.01A	10	10/19	0.2	0.10	2
142.202.02A	25	14/23	0.5	0.25	2
142.202.03A	50	19/26	1.0	0.50	2
142.202.04A	100	24/29	1.0	0.50	2
142.202.05A	250	29/32	2.0	1.00	2
142.202.06A	500	34/35	5.0	2.50	2
142.202.07A	1000	45/40	10.0	5.00	1
142.202.08A	2000	45/40	20.0	10.00	1

3.3  
BORO  
GLASS



A  
CLASS



ENQUIRE

\* Non DIN



## Mixing Cylinder with Hexagonal base, Class-B

- Complies with DIN EN ISO 4788 standards
- With PE stopper

Item Code	Capacity (ml)	N/S	Sub Division (ml)	Tolerance ( $\pm$ ml)	Pack of
142.222.00A*	5	10/19	0.1	0.1	2
142.222.01A	10	10/19	0.2	0.2	2
142.222.02A	25	14/23	0.5	0.5	2
142.222.03A	50	19/26	1.0	1.0	2
142.222.04A	100	24/29	1.0	1.0	2
142.222.05A	250	29/32	2.0	2.0	2
142.222.06A	500	34/35	5.0	5.0	2
142.222.07A	1000	45/40	10.0	10.0	1
142.222.08A	2000	45/40	20.0	20.0	1

3.3  
BORO  
GLASS



B  
CLASS



ENQUIRE

\* Non DIN