

# Flask



## With Hollow Glass Stopper

Item Code	Capacity (ml)	Socket Size	Height (mm)	Approx Dia (mm)	Pack of
076.202.01A	100	24/29	64	100	2
076.202.01	100	29/32	64	100	2
076.202.02	250	29/32	85	140	2
076.202.03A	500	24/29	105	175	2
076.202.03	500	29/32	105	175	2
076.202.04	1000	29/32	131	220	2
076.202.05	2000	29/32	166	270	1

3.3  
BORO  
GLASS



ENQUIRE



## Flask, Conical (Erlenmeyer) with Screw Cap

- Made from heat resistant ASTM E-438 Type 1, Boro 3.3 glass
- With autoclavable PP screw cap and PTFE Liner

Item Code	Capacity (ml)	Thread Size	Dia (mm)	Pack of
075.202.01	50	GL 25	51	2
075.202.02	100	GL 25	64	2
075.202.03	150	GL 25	72	2
075.207.04	250	GL 45	85	2
075.207.05	500	GL 45	105	2
075.207.06	1000	GL 45	131	2
075.207.07	2000	GL 45	166	2
075.207.10	5000	GL 45	220	2

3.3  
BORO  
GLASS



ENQUIRE