

# Funnels



**Funnel, Dropping, Gilson With Glass Stopcock**  
• With PE Stopper

Item Code	Capacity (ml)	N/S	Bore (mm)	Stopcock Type	Pack of
159.202.01	50	19/26	2.5	Solid Glass	1
159.202.02	100	19/26	2.5	Solid Glass	1
159.202.03	250	29/32	4	Solid Glass	1
159.202.04	500	29/32	4	Solid Glass	1
159.202.05	1000	29/32	6	Solid Glass	1
159.202.06	2000	29/32	6	Solid Glass	1

3.3  
BORO  
GLASS

Chemical  
Resistant

WHITE  
PRINTING

GLASS

PE



ENQUIRE



**Funnel, Dropping, Gilson With PTFE Stopcock**  
• With PE Stopper

Item Code	Capacity (ml)	N/S	Bore Size	Stopcock Type	Pack of
160.202.01	50	19/26	2.5	PTFE	1
160.202.02	100	19/26	2.5	PTFE	1
160.202.03	250	29/32	4	PTFE	1
160.202.04	500	29/32	4	PTFE	1
160.202.05	1000	29/32	6	PTFE	1
160.202.06	2000	29/32	6	PTFE	1

3.3  
BORO  
GLASS

Chemical  
Resistant

WHITE  
PRINTING

PE



ENQUIRE