

Funnels



Funnel, Separating, Squibb shape With PTFE Stopcock

- Complies with ISO 4800 Standard
- Without Graduation
- With PE Stopper

Item Code	Capacity (ml)	N/S	Stopcock Type	Bore Size	Pack of
149.202.02	50	19/26	PTFE	2.5	1
149.202.02A*	60	14/23	PTFE	2.5	1
149.202.03	100	19/26	PTFE	2.5	1
149.202.03A*	125	19/26	PTFE	2.5	1
149.202.04	250	29/32	PTFE	4	1
149.202.05	500	29/32	PTFE	4	1
149.202.06	1000	29/32	PTFE	6	1
149.202.07	2000	29/32	PTFE	6	1

3.3
BORO
GLASS

Chemical
Resistant

WHITE
PRINTING

PTFE

PTFE

PE



ENQUIRE

* Non DIN



Funnel, Separating, Squibb shape With PTFE Stopcock

- Complies with ISO 4800 Standard
- With Graduation
- With PE Stopper

Item Code	Capacity (ml)	N/S	Stopcock Type	Bore Size	Pack of
149.209.02	50	19/26	PTFE	2.5	1
149.209.03	100	19/26	PTFE	2.5	1
149.209.04	250	29/32	PTFE	4	1
149.209.05	500	29/32	PTFE	4	1
149.209.06	1000	29/32	PTFE	6	1
149.209.07	2000	29/32	PTFE	6	1

3.3
BORO
GLASS

Chemical
Resistant

WHITE
PRINTING

PTFE

PTFE

PE



ENQUIRE