



Plastic Pasteur pipettes

Code DLB010, DLB011, DLB022, DLB023, DLB024, DLB025, DLB026, DLB027, DLB028, DLB029, DLB030, DLB031, DLB032, DLB033

- Made of transparent low-density polyethylene (LDPE)
- DNase/RNase, endotoxin and pyrogen free
- Good chemical stability
- Embossed graduation (except DLB032 and DLB033)
- Accurate and with good repeatability
- Individually wrapped sterile pipettes. Sterilised with ethylene oxide.



Code	Packaging
DLB010	500 pcs/box
DLB011	500 pcs/box
DLB022	1000 pcs./bag
DLB023	1000 pcs./bag
DLB024	500 pcs./bag
DLB025	500 pcs./bag
DLB026	500 pcs./bag
DLB027	500 pcs./bag
DLB028	100 pcs./bag
DLB029	100 pcs./bag
DLB030	100 pcs./bag
DLB031	100 pcs./bag
DLB032	500 pcs/bag
DLB033	500 pcs/bag

Technical specifications

Code	Nominal volume (mL)	Total volume (mL)	Drop volume (μL)	Total length (mm)	Sterile
DLB010	1	5	45	145	No
DLB011	3	7	48	155	No
DLB022	0,2	0,9	20	65	No
DLB023	0,5	1,8	35	126	No
DLB024	1	3	42	145	No
DLB025	2	5	48	155	No
DLB026	3	6	48	155	No
DLB027	3 (long)	8,2	48	180	No
DLB028	5	12	40	210	No
DLB029	10	23	55	300	No
DLB030	1	5	45	145	Yes
DLB031	3	7	48	155	Yes
DLB032	1	5,8	15	149	No
DLB033	1	5,8	15	149	Yes