



Scalpels

- Scalpels are very sharp knives used in surgical operations, anatomic dissections and autopsies, being very useful when a precise and thin cut is required.
- The scalpels present two components, blade and handle, that are available in several sizes and shapes; the use of each model will depend on the type and accessibility of the tissues to be cut.



Disposable scalpels

From FMB004 to FMB011

- They are supplied in boxes of 10 units. Each scalpel is complete with blade made of carbon steel and plastic handle, totally sterile and individually packaged.
- With CE 0434 certificate.

Code	Number	Length
FMB004	10	135 mm
FMB005	11	135 mm
FMB006	15	135 mm
FMB007	20	160 mm
FMB008	21	160 mm
FMB009	22	160 mm
FMB010	23	160 mm
FMB011	24	160 mm



Scalpel blades

From FMB012 to FMB019, FMB026

- Made of carbon steel and extremely sharp, they require a minimum pressure to make a precise cut.
- Scalpel blades differ in their size and shape, and according to these features each model is represented by a universal number.
- They are supplied in boxes of 100 units totally sterile and individually packaged.
- With CE 0434 certificate.

Code	Number	For handle number
FMB012	Nº 10	3
FMB013	Nº 11	3
FMB026	Nº 13	3
FMB014	Nº 15	3
FMB015	Nº 20	4
FMB016	Nº 21	4
FMB017	Nº 22	4
FMB018	Nº 23	4
FMB019	Nº 24	4



Scalpel handles

- To be used with scalpel blades. Nahita handles provide a firm holding of the blade, without risk of release.
- Made of stainless steel, each model is available in two sizes.



N°3

N° 3

FMB020, FMB021

Code	Number	Length
FMB020	N° 3 long	200 mm
FMB021	N° 3	125 mm



N°4

N° 4

FMB022, FMB023

Code	Number	Length
FMB022	N° 4 long	220 mm
FMB023	N° 4	140 mm