

CHROMATOGRAPHY FILTER PAPER

Code PKC001, PKC003, PKC004, PKC006, PKC007, PKC008, PKC009, PKC010

Chromatography filter paper designed for analytical applications requiring high reproducibility and precise control of the separation process. Available in different dimensions, the paper presents the following properties:

- Basis weight: 76–84 g/m²
- Filtration time: 35–70 s
- Wet bursting strength: ≥ 50 kPa
- Dry bursting strength: ≥ 120 kPa
- Whiteness: ≥ 82%
- pH: 6–8



Code	Dimensions	Uds/box
PKC001	100x300 mm	100
PKC006	140x60 mm	100
PKC003	200x200 mm	100
PKC004	250x250 mm	100
PKC007	460x570 mm	100
PKC008	580x600 mm	100
PKC009	580x600 mm	25
PKC010	285x250 mm	25