


# Flask



## Flask, Conical (Erlenmeyer), Heavy Duty with Narrow Mouth

- Made from heat resistant ASTM E-438 Type 1, Boro 3.3 glass
- Complies with ISO 1773 standard
- All the flasks are printed in durable white enamel
- Has a reinforced thick rim to reduce chipping & breakage

Item Code	Capacity (ml)	Body Dia (mm)	Neck Dia (mm)	Height (mm)	Pack of
231.237.01	25	42	22	75	10
231.237.02	50	51	22	90	10
231.237.03	100	64	22	105	10
231.237.03A*	150	74	28	118	10
231.237.05	250	85	34	145	10
231.237.07	500	105	34	180	10
231.237.08	1000	131	42	220	10
231.237.09	2000	166	50	280	6
231.237.10	3000	187	50	310	1
231.237.12	5000	220	50	365	1

3.3  
BORO  
GLASS

WHITE  
PRINTING


Chemical  
Resistant

HEAVY  
WALL



ENQUIRE

\* Non DIN



## Flask, Conical (Erlenmeyer) Wide Mouth

- Made from heat resistant ASTM E-438 Type 1, Boro 3.3 glass
- Complies with ISO 24450 standard
- All the flasks are printed in durable white enamel
- Has a reinforced thick rim to reduce chipping & breakage

Item Code	Capacity (ml)	Dia (mm)	Dia1 (mm)	Height (mm)	Pack of
232.202.01*	25	42	31	70	10
232.202.02	50	51	34	85	10
232.202.03	100	64	34	105	10
232.202.05	250	85	50	140	10
232.202.07	500	105	50	175	10
232.202.08	1000	131	50	220	10
232.202.09*	2000	166	72	276	6

3.3  
BORO  
GLASS

WHITE  
PRINTING

Chemical  
Resistant



ENQUIRE

\* Non DIN