



ROUND ULTRASONIC BATHS WITH HEATING

Code ZCC031, ZCC032, ZCC033

- AISI 304 stainless steel tank: corrosion-resistant and highly durable material, suitable for industrial use
- Industrial-grade integrated circuit ensures stable and reliable operation during prolonged use
- Temperature range: 0-80°C
- Adjustable operating time: 1-99 minutes
- Wax removal and degreasing functions: suitable for removing oils, waxes and stubborn residues
- Operating frequency:
 - 28 kHz for general cleaning or pre-washing
 - 40 kHz for intensive or precision cleaning
- Automatic standby and sleep modes, activated with a single key
- High-precision, long-life digital controller facilitates accurate adjustment of operating parameters and easy maintenance



Technical Specifications:

Code	ZCC031	ZCC032	ZCC033
Tank dimensions	Ø300x220 mm	Ø350x240 mm	Ø450x280 mm
Overall dimensions	410x350x520 mm	460x400x540 mm	560x500x580 mm
Capacity	15 L	22 L	44 L
Ultrasonic power	360 W	480 W	900 W
Heating power	400 W	600 W	1000 W
Frequency	28/40 kHz		
Time	1-99 min		
Temperature	0-80 °C		
Voltage	AC 220V		