

Beakers



Beaker, Griffin/Low Form With Double Capacity Scale, Heavy Wall

- Complies with ASTM E960, TYPE-2, BORO 3.3 GLASS
- Uniform wall thickness with reinforced rim and spout
- Printed in matt finish white colour enamel and double scale
- Highly resistant to chemicals and High mechanical strength
- Not recommended to use on hotplates

| Item Code | Capacity (ml) | Graduation Interval (ml) | Graduation Range (ml) | Pack of |
|-------------|---------------|--------------------------|-----------------------|---------|
| 229.235.07 | 250 | 25 | 25-200 | 12 |
| 229.235.08 | 400 | 25 | 25-325 | 12 |
| 229.235.09 | 600 | 50 | 50-500 | 6 |
| 229.235.10 | 1000 | 50 | 50-1000 | 6 |
| 229.235.11* | 2000 | 100 | 200-1800 | 4 |
| 229.235.14 | 4000 | 250 | 250-3500 | 1 |

3.3
BORO
GLASS

Chemical
Resistant

WHITE
PRINTING

8 I
LOW FORM

2x METRIC

HEAVY
WALL



ENQUIRE

*Non ASTM Standards



Beaker, Berzelius/Tall Form With Spout

- Complies with DIN 12331 and ISO 3819 standards
- Clear metric Graduation for convenient reading
- Highly resistant to chemicals and thermal shock
- Large labeling field for easy marking
- Uniform wall thickness and high mechanical strength

| Item Code | Capacity (ml) | Dia (mm) | Height (mm) | Pack of |
|--------------|---------------|----------|-------------|---------|
| 230.202.01 | 25 | 32 | 60 | 10 |
| 230.202.02 | 50 | 38 | 70 | 10 |
| 230.202.03 | 100 | 48 | 80 | 10 |
| 230.202.04 | 150 | 54 | 95 | 10 |
| 230.202.05 | 250 | 60 | 120 | 10 |
| 230.202.06 | 400 | 70 | 130 | 10 |
| 230.202.07A* | 500 | 75 | 142 | 8 |
| 230.202.07 | 600 | 80 | 150 | 8 |
| 230.202.08 | 1000 | 95 | 180 | 6 |

3.3
BORO
GLASS

Chemical
Resistant

WHITE
PRINTING

8 I
TALL FORM



ENQUIRE

*Non DIN