



Distillation flasks

Code PBG001, PBG002, PBG003, PBG004

- As the name suggests, distillation flasks are used in reactions that require the evaporation of liquids and subsequent condensation.
- For this purpose, Endo Glassware distillation flasks are made of borosilicate glass, highly resistant to thermal shock, and feature a spherical base for uniform heating of the solutions to be distilled.
- They also have a long, narrow neck and a side arm for connection to a condenser, allowing the condensation of the evaporated liquid.



Referencia	Capacity (mL)	Ø Neck (mm)	Height (mm)	Ø Arm (mm)	Arm Length (mm)
PBG001	100	15	174	7	100
PBG002	250	19	216	8	120
PBG003	500	24	275	8	150
PBG004	1000	36	350	9	250