



# AUXILAB BENCHTOP AUTOCLAVES

Code ZCK016, ZCK017, ZCK018, ZCK019

- Small, fully automatic class B autoclaves
- Easy operation and maintenance
- LCD screen
- Interface available in Spanish and English
- Six predefined sterilisation programmes for solid, porous and hollow materials
- Special predefined sterilisation programme for prion inactivation
- User-defined programme
- Programmes for vacuum, BD (Bowie-Dick) and Helix tests
- Cleaning programme
- Option to set a user password
- With printer and USB output
- Electrical and mechanical door locking system
- Electrical protection against short circuits



Code	ZCK016	ZCK017	ZCK018	ZCK019
Model	AXL8L	AXL12L	AXL18L	AXL22L
Capacity	8 L	12 L	18 L	22 L
Sterilisation temperature	105-134°C			
Sterilisation time	4-60 min			
Maximum working pressure	0,21-0,23 MPa			
Vacuum cycles	1-10			
Drying time	1-25 min			
Chamber dimensions	Ø170x320 mm	Ø200x360 mm	Ø250x350 mm	Ø250x450 mm
Chamber material	304 stainless steel			
Equipment dimensions (LxWxH)	640 x 440 x 365 mm		670 x 500 x 435 mm	
Net weight	45 kg	48 kg	53 kg	55 kg
Power supply	220V, 50Hz			
Power	1,95 kW			