



CERTIFICATE OF CONFORMITY

Issuer's Name / Producer : GLASSCO LABORATORY EQUIPMENTS PVT. LTD.
Issuer's Address / Producer : Manglai, P.O. Khudda Kalan, Ambala Cantt - 133001
Haryana, India
Object of Declaration : Tubes,
Material : Borosilicate Glass
Purpose of Use : Tubes used by chemists to hold, mix, or heat small quantities of solid or liquid chemicals, especially for qualitative experiments and assays.

Specification of Tubes specified in Table-1, as per Glassco / ISO / or ASTM Standard

Catalogue No.	Product Description	Bottom Type	O.D. x Length	Pack Qty	
082.202.01	Round Bottom Tube with screw cap	Round	15 x 75	100	
082.202.02	Round Bottom Tube with screw cap	Round	15 x 125	100	
082.202.03	Round Bottom Tube with screw cap	Round	15 x 150	100	
082.202.04	Round Bottom Tube with screw cap	Round	18 x 150	100	
082.202.05A	Round Bottom Tube with screw cap	Round	25 x 95	100	
082.202.05B	Round Bottom Tube with screw cap	Round	25 x 50	100	
082.202.06	Round Bottom Tube with screw cap	Round	25 x 150	100	

The object of declaration described above is in conformity with the requirements of the following documents:

- **DIN ISO 3585 Borosilicate Glass 3.3**

Material: The above are made from DIN ISO 3585:1998, heat resistant Boro 3.3 Glass that meet ASTM E-438, TYPE-1, CLASS- A Standard.

Compliance: These Tubes comply in every respect characteristics, measures and accuracy are in full conformity With the provisions of the standard ISO / ASTM / or Glassco Specification.

Additional Information:

Glassco has installed the quality control system as per ISO 9001:2008 and products delivered to the market are in full conformity with the technical Specification as per standard and with all fundamental requirements to such products.

Certificate No. IPK 0039.69-01 issued by intertek certification limited, A UKAS Management System Accredited Body



This document has been produced electronically & is valid without signature