

CERTIFICATE OF CONFORMITY

Descripción | Description | Description: **FUNNEL WITH HANDLE, STAINLESS STEEL 316**

Referencia | Code | Référence: **PNN004, PNN005, PNN006, PNN007, PNN008**

Marca | Brand | Marque: **Endo Metalware**

Material | Material | Matériel: **AISI-316**

Auxilab declares that the product described above has been manufactured and controlled under a Quality Management System certified according to ISO 9001:2015.

The material has been tested by Optical Emission Spectrometry (OES) to verify its chemical composition according to the applicable material specification. The test results are expressed as percentage by mass and are shown below.

Chemical Composition	
Element	Result (%)
Cr	16.129
Mn	1.318
S	0.008
Mo	2.053
C	0.031
P	0.011
Si	0.453
Ni	10.272

Fecha | Date | Date: 22/04/2026

Sello | Stamp | Timbre:

 **AUXILAB, S. L.**
N.I.F. ESB31072210
MATERIAL PARA LABORATORIO
Polig. Morea Norte, 8 - Tel. 34-948 310 513
31191 BERRIAIN - NAVARRA (ESPAÑA)