



Sterile handles

Code **BRG001, BRG002, BRG003, BRG024**

- Manufactured from medical grade polystyrene.
- Highly resistant, suitable for inoculation and scratching.
- Handle diameter ensures inoculation accuracy.
- Ethylene oxide sterilised.
- In bags of 20 or individually packaged.



Technical specifications

| Code | BRG001 | BRG002 | BRG003 | BRG024 |
|----------------------|--------|--------|--------|--------|
| Volume (µl) | 1 | 10 | 10 | 1 |
| Format (units/bag) | 20 | 1 | 20 | 1 |
| Packaging (bags/lot) | 25 | 200 | 25 | 200 |
| Lot | 500 | 100 | 500 | 100 |

Sterile handles for extensions

Code **BRG022, BRG023**

- Made of ABS
- L-shape
- Dimensions (mm): 144x33xØ4
- Ethylene oxide sterilised.
- In bags of 10 pieces or individual packaging.



Technical specifications

| Code | BRG022 | BRG023 |
|----------------------|--------|--------|
| Format (units/bag) | 10 | 1 |
| Packaging (bags/lot) | 25 | 100 |
| Lot | 250 | 100 |

Kolle Handle

Code **BRG009, BRG010**

- Used for surface or deep seeding of bacterial cultures.
- Each handle is supplied with two flame-sterilisable Ni-Cr seeding handles which are easily inserted into the end of the handle.



Technical specifications

| Code | BRG009 | BRG010 |
|-------------|---------|----------|
| Length (mm) | 200 | 290 |
| Handle | Plastic | Metallic |



Kolle Handle Accessories

Code **BRG011, BRG012, BRG013**

- Made of Nickel-Chrome



Technical specifications

| Code | BRG011 | BRG012 | BRG013 |
|-------------------|-----------------|-------------------|------------|
| Description | Ø 3 mm ring end | Needle-shaped end | Lancet tip |
| Total length (mm) | 43 | 48 | 59 |
| Lot | 50 | 50 | 59 |

L-handles for extensions

Code **BRG017, BRG018**

- Made of stainless steel
- Reusable
- Autoclavable and flame sterilisable



Technical specifications

| Code | BRG017 | BRG018 |
|-------------------|--------|--------|
| Head (mm) | 30 | 138 |
| Total length (mm) | 150 | 150 |

Triangular handles for extensions

Code **BRG019, BRG020, BRG021**

- Made of stainless steel
- Reusable
- Autoclavable and flame sterilisable



Technical specifications

| Code | BRG019 | BRG020 | BRG021 |
|-------------------|--------|--------|--------|
| Head (mm) | 16 | 25 | 30 |
| Total length (mm) | 200 | 200 | 200 |