




# CERTIFICATE OF CONFORMITY

Issuer's Name / Producer : GLASSCO LABORATORY EQUIPMENTS PVT. LTD.  
Issuer's Address / Producer : Manglai, P.O. Khudda Kalan, Ambala Cantt - 133001  
Haryana, India  
Object of Declaration : Measuring Cylinder  
Material : Borosilicate Glass  
Purpose of Use : Measuring Cylinder is used in a laboratory for measuring exact quantities of a liquid. It is also called a graduated cylinder, as it is marked with precise measurements.

## Specification of cylinder specified in Table-1, EN ISO 4788 : 2005 (D)

Catalogue No.	Product Description	Capacity mL.	Pack Qty	
139.PK1.01A	Measuring Cylinder, Hexagonal Base Batch Certified	10	1	
139.PK1.02A	Measuring Cylinder, Hexagonal Base Batch Certified	25	1	
139.PK1.03A	Measuring Cylinder, Hexagonal Base Batch Certified	50	1	
139.PK1.04A	Measuring Cylinder, Hexagonal Base Batch Certified	100	1	
139.PK1.05A	Measuring Cylinder, Hexagonal Base Batch Certified	250	1	
139.PK1.06A	Measuring Cylinder, Hexagonal Base Batch Certified	500	1	

## The object of declaration described above is in conformity with the requirements of the following documents:

- DIN ISO 3585 Borosilicate Glass 3.3
- EN ISO 4788 Laboratory Glassware - Measuring Cylinder

**Material:** The above are made from DIN ISO 3585:1998, heat resistant Boro 3.3 Glass that meet ASTM E-438, TYPE-1, CLASS- A Standard.

**Compliance:** These cylinders comply in every respect characteristics, measures and accuracy are in full conformity with the provisions of the standard EN ISO 4788: 2005 (D)

## Additional Information:

Glassco has installed the quality control system as per ISO 9001 and products delivered to the market are in full conformity with the technical Specification as per standard and with all fundamental requirements to such products.

Certificate No. IPK 0039.69-01 issued by intertek certification limited, A UKAS Management System Accredited Body



This document has been produced electronically & is valid without signature