




# CERTIFICATE OF CONFORMITY

Issuer's Name / Producer : GLASSCO LABORATORY EQUIPMENTS PVT. LTD.  
Issuer's Address / Producer : Manglai, P.O. Khudda Kalan, Ambala Cantt - 133004  
Haryana, India  
Object of Declaration : Cones with Screw thread  
Material : Borosilicate Glass  
Purpose of Use : This threaded joints provide a vacuum tight seal without grease  
frozen or sized joints, can be taken apart using the loosening ring.

## Specification of Cone specified in Table-1, as per Glassco specification

Catalogue No.	Product Description	Socket size	Pack Qty.
002.SC.04	CONES WITH SCREW THREAD	14/23	10
002.SC.05	CONES WITH SCREW THREAD	19/26	10
002.SC.06	CONES WITH SCREW THREAD	24/29	10
002.SC.07	CONES WITH SCREW THREAD	29/32	10
002.SC.08	CONES WITH SCREW THREAD	34/35	10
002.SC.09A	CONES WITH SCREW THREAD	45/40	10



## The object of declaration described above is in conformity with the requirements of the following documents:

- **DIN ISO 3585** **Borosilicate Glass 3.3**

**Material:** The above are made from DIN ISO 3585:1998, heat resistant Boro 3.3 Glass that meet ASTM E-438, TYPE-1, CLASS- A Standard.

## Additional Information:

Glassco has installed the quality control system as per ISO 9001 and products delivered to the market are in full conformity with the technical Specification as per standard and with all fundamental requirements to such products.

Certificate No. IPK 0039.06 issued by intertek certification limited, A UKAS Management System Accredited Body



This document has been produced electronically & is valid without signature