



## Vacuum sealer "Sealcom-V"

OPTIMUM QUALITY SEAL AND PRESERVATION.  
SUITABLE FOR SEALING SOLIDS AND LIQUIDS.

### APPLICATIONS

Vacuum seal for any type of thermoplastic material: polyethylene, polypropylene, polyester, aluminium etc. Ideal for sealing solids, liquids, powders, granules, seeds, fruit, food, medical apparatus, forensics, chemical products, sludge samples, mud, cores, etc. Prevents the oxidation, mould and decomposition. Prolongs the quality, durability and freshness of items to be stored.

### FEATURES

Bench top. Stainless steel AISI 304 robust construction with transparent lid.  
Easy to use and program. After closing the door the complete process is automatic.  
Operates using a oil vacuum pump (type AV-30 or N32).

### CONTROL PANEL

Vacuum gauge with scale 0 to -0.1 Mpa.  
Sealer temperature adjustment, high, medium or low.  
Push button emergency stop.  
LED indicator of operation.  
Digital selection and display of the following parameters:  
- Vacuum time 0-99 seconds.  
- Seal time 0 to 9.9 seconds  
- Discharge time 0 to 9.9 seconds



MODEL	Part No.	Seal Width mm	Max. seal length mm	Pump rate m <sup>3</sup> /h	Vacuum limit Mpa	Height / Width / Depth (usable) cm	Height / Width / Depth o (exterior) cm	Power W	Weight Kg
SEALCOM-V	4100050	10	260	10	-0.1	5 28 39	36 35 52	370	35

Supplied with a maintenance kit that includes 2 heaters, 3 PTFE heater covers, 1 lid seal and a 250 ml bottle of oil.

**ADDITIONAL SPARES** Maintenance kit for Sealcom-V: 2 heaters, 3 PTFE heater covers, 1 lid seal. Part No. 4100051  
Vacuum pump oil: 1 litre bottle (Type AV-30), Part No. 4100052



## Bag sealer "S-400-B"

SUITABLE FOR POLYPROPYLENE-POLYESTER STERILIZING BAGS WITH A PAPER BASE.

### FEATURES

Metal epoxy covered case.  
Pressure sealer. Manual cut.

### CONTROL PANELS

Illuminated mains switch.  
Lamp "ON" when the thermostat has reached the pre-set working temperature.

### MODEL

Part No.	Seal Width mm	Maximum seal length mm	Height / Width / Depth (exterior) cm	Power W	Weight Kg
4100047	12	300	12 35 21	100	4,8

### SEALING MACHINE ACCESSORIES

#### Rolls of paper for sterilisation "Polisteril"

Paper base polypropylene and polyester bags showing differing ink prints of sterilized items.

Part No.	Bag width cm	Bag length m
4100075	7,5	100
4100100	10	100
4100200	20	100



SEE A FULL RANGE OF BAG SEALERS FOR STERILIZATION AND ACCESSORIES ON PAGE 68 AND 69

Automatic bag sealer  
"Seal-Com 2". Part No: 4100044



Bag sealer  
"Seal-1000". Part No: 4101712

