



Laboratory Jack

Designed to lift and maintain all types of equipment and apparatus such as flasks, flasks, baths, stirrers, heating plates, etc.

Laboratory jack aluminium

Code FCD001, FCD002, FCD003



- Made of aluminium
- With scissor-type support
- With manual control and non-slip base
- They have 4 non-slip rubber feet to ensure good stability and avoid sudden movements.

Code	Platform (mm)	Height (mm)	Max. load (Kg)
FCD001	100x100	40-140	14
FCD002	150x150	50-270	27
FCD003	200x200	50-320	38

Laboratory jack stainless steel

Code FCD004, FCD005, FCD006, FCD007, FCD008



- Made of AISI 304 stainless steel
- With scissor-type support
- With manual control
- They have 4 non-slip rubber feet to ensure good stability and avoid sudden movements.

Code	Platform (mm)	Height (mm)	Max. load (Kg)
FCD004	100x100	55-140	45
FCD005	150x150	65-250	60
FCD006	200x200	65-250	60
FCD007	250x250	70-385	75
FCD008	300x300	70-385	75