



## Double bent connector, ENDO GLASSWARE

Code PBQ021, PBQ022, PBQ034, PBQ035

- Made of borosilicate glass



Code	Ø (mm)	Length (mm)	Frosted
PBQ021	7	160	-
PBQ022	12	250	-
PBQ034	7	160	14/23 (M/M)
PBQ035	25	280	29/32 (M/M)

## Thermometer adapter with screw cap, ENDO GLASSWARE

Code PBQ026, PBQ027, PBQ028

- Made of borosilicate glass
- With threaded head for 5 to 7 mm thermometers



Code	Thread	Frosted	Angle (°)
PBQ026	13	14/23 (M)	-
PBQ027	13	14/23 (M/M)	75
PBQ028	18	29/32 (M/M)	75

## Connecting adapter with hose, ENDO GLASSWARE

Code PBQ029

- Made of borosilicate glass
- For laboratory glass assemblies



Code	Frosted	Ø hose connection (mm)
PBQ029	14/23 (M/H)	10

## Bent distillation adapter joint, ENDO GLASSWARE

Code PBQ031

- Made of borosilicate glass
- Short stem



Code	Frosted	Stem length (mm)	Angle (°)	Ø ext. outlet pipe (mm)
PBQ031	29/32 (H)	100	105	12

## Thermometer Claisen adapter, ENDO GLASSWARE

Code PBQ032

- Made of borosilicate glass
- With threaded head for 5 to 7 mm thermometers




Code	Thread	Frosted
PBQ032	13	14/23 (H/M/M)

## Receiver adapter for 2 flasks, ENDO GLASSWARE

Code PBQ033

- Made of borosilicate glass
- For two flasks



Code	Ø hose connection (mm)	Frosted
PBQ033	10	14/23 (H/M)