



## Pipeclay triangles

**Code FFG001, FFG002, FFG003, FFG004**

Nahita pipeclay triangles are made up of three refractory sillimanitebased clay tubes which are joined together in the shape of a triangle by a galvanized metal wire.

Given their high resistance to extreme temperatures without cracking or deforming, they are used for supporting crucibles on a tripod during the calcination process.



Code	Side length	Tube Ø
FFG001	30 mm	9 mm
FFG002	50 mm	8 mm
FFG003	65 mm	8 mm
FFG004	75 mm	8 mm