



AXAVET OLV-5 Oxygen Concentrator

Code **ZMC017**

The AXAVET 5 oxygen concentrator has been designed to meet the needs of veterinary centres requiring oxygen supply for anaesthesia, oxygen therapy and nebulisation. Its output pressure is higher than most other concentrators on the market, allowing it to efficiently deliver oxygen to all brands and models of anaesthesia machines.



Features

- Fully automatic and easy operation
- Compact and silent
- Requires very little space
- With supply cut-off alarm
- Includes nebuliser cup
- Easy maintenance

Technical specifications

Code	ZMC017
Oxygen purity	$\geq 93\% \pm 3\%$
Output pressure	100-115 kPa
Flow rate range	0-5 L/min
Noise level	< 45 dB
Power consumption	500 W
Power supply	100-240 V AC, 50/60 Hz
Dimensions (LxWxH)	490x270x560 mm
Weight	22 kg