

Beakers



Beaker, Griffin/Low Form With Double Capacity Scale, Heavy Wall

- Complies with ASTM E960, TYPE-2, BORO 3.3 GLASS
- Uniform wall thickness with reinforced rim and spout
- Printed in matt finish white colour enamel and double scale
- Highly resistant to chemicals and High mechanical strength
- Not recommended to use on hotplates

Item Code	Capacity (ml)	Graduation Interval (ml)	Graduation Range (ml)	Pack of
229.235.07	250	25	25-200	12
229.235.08	400	25	25-325	12
229.235.09	600	50	50-500	6
229.235.10	1000	50	50-1000	6
229.235.11*	2000	100	200-1800	4
229.235.14	4000	250	250-3500	1

3.3
BORO
GLASS

Chemical
Resistant

WHITE
PRINTING

8 I
LOW FORM

2x METRIC

HEAVY
WALL



ENQUIRE

*Non ASTM Standards



Beaker, Berzelius/Tall Form With Spout

- Complies with DIN 12331 and ISO 3819 standards
- Clear metric Graduation for convenient reading
- Highly resistant to chemicals and thermal shock
- Large labeling field for easy marking
- Uniform wall thickness and high mechanical strength

Item Code	Capacity (ml)	Dia (mm)	Height (mm)	Pack of
230.202.01	25	32	60	10
230.202.02	50	38	70	10
230.202.03	100	48	80	10
230.202.04	150	54	95	10
230.202.05	250	60	120	10
230.202.06	400	70	130	10
230.202.07A*	500	75	142	8
230.202.07	600	80	150	8
230.202.08	1000	95	180	6

3.3
BORO
GLASS

Chemical
Resistant

WHITE
PRINTING

8 I
TALL FORM



ENQUIRE

*Non DIN