



## Cell culture flasks

Code **CGR001, CGR002, CGR003, CGR004, CGR005, CGR006**

- Manufactured from 100% virgin polystyrene with high transparency
- Surface treated (TC) for better cell adhesion
- Electron beam sterilised
- Non-pyrogenic and DNase/RNase free
- Available with standard and vent caps (with 0.22  $\mu\text{m}$  hydrophobic filters)
- With growth areas of 25, 75 and 175  $\text{cm}^2$ , and volumes of 50, 250 and 600 mL, respectively
- Graduations on both sides of the flask
- Stackable
- Lot number indicated on the flask for traceability



Reference	Packaging
CGR001	5 pcs/bag, 40 bags/carton
CGR002	5 pcs/bag, 40 bags/carton
CGR003	5 pcs/bag, 20 bags/carton
CGR004	5 pcs/bag, 20 bags/carton
CGR005	5 pcs/bag, 10 bags/carton
CGR006	5 pcs/bag, 10 bags/carton

## Technical specifications

Reference	CGR001	CGR002	CGR003	CGR004	CGR005	CGR006
Growth area ( $\text{cm}^2$ )	25	25	75	75	175	175
Volume (mL)	50	50	250	250	600	600
Cap type	Vent	Standard	Vent	Standard	Vent	Standard