

Adapter



Adapter, Receiver Plain Bend with Vent

- With 105° bend drip tip
- Sidearm vent provides an escape to gases during experiments and avoids closed system problems
- Used for distribution of liquid from the condenser into collection vessel

Item Code	Socket Size	Cone Size	Pack of
021.202.01	14/23	14/23	1
021.202.02	19/26	19/26	1
021.202.03	19/26	24/29	1
021.202.04	24/29	24/29	1
021.202.05	29/32	29/32	1

3.3
BORO
GLASS



1 x
SOCKET

1 x
CONE



ENQUIRE



Adapter, Receiver Plain Bend without Vent

- With 105° bend drip tip
- Without Vent
- Used for distribution of liquid from the condenser into collection vessel

Item Code	Socket Size	Cone Size	Pack of
025.202.00	14/23	14/23	1
025.201.00A	14/23	19/26	1
025.202.01	19/26	19/26	1
025.202.02	19/26	24/29	1
025.202.03	24/29	24/29	1
025.202.04	29/32	29/32	1

ASTM Sizes

Item Code	Socket Size	Cone Size	Pack of
025.470.06	24/40	24/40	1

3.3
BORO
GLASS



1 x
SOCKET

1 x
CONE



ENQUIRE



ENQUIRE



Adapter, Distillation Angled 105°

- With drip tip and delivery stem inside the cone
- Used for distribution of liquid from the condenser into collection vessel

Item Code	Socket Size	Cone Size	Pack of
025.475.01	14/20	14/20	1

3.3
BORO
GLASS



1 x
SOCKET

1 x
CONE



ENQUIRE