



## Water and oil baths “Termobat”

ADJUSTABLE TEMPERATURE FROM AMBIENT+5 °C TO 110 °C. STABILITY ±2 °C TO 200 °C. STABILITY ±3 °C.

### SAFETY:

CONFORMS TO THE DIN STANDARD 12877 BY INCORPORATING A MANUALLY RESET OVER TEMPERATURE SAFETY CUT-OUT.

### APPLICATIONS

To be used with a rotary evaporator for flasks up to 10 litres or as a stand alone bain marie.

### FEATURES

Tank interior made of AISI 304 stainless steel, with an epoxy coated base. The heater element is embedded in a protective INCOLOY stainless steel shield.

### CONTROL PANEL

Main power switch, with luminous indicator lamp. A hydraulic temperature thermostat regulator is synchronised with an indicator lamp.

### MODELS

Part No.	Maximum Temperature °C	Capacity litres	Ø (usable) cm	Height cm	Ø (exterior) cm	Height cm	Power W	Weight Kg
3000587	110	11	29	17	29	25	2100	5.5
3000588	200	11	29	17	29	25	2100	5.5

### Large capacity



### ACCESSORIES MADE OF AISI 304 STAINLESS STEEL

Reduction ring cover lid of 210, 180, 145, 115 and 80 mm Ø. With a location hole for a thermometer. Part No. 7000589



Heater element cover shelf. Part No. 6000055



Thermometer. Part No.

1095100 up to 100 °C.  
1095250 up to 250 °C.



## Sand baths “Combiplac-Sand”

ADJUSTABLE PLATE TEMPERATURE UP TO 400 °C.



### FEATURES

Electronic energy heater controller. The steel hotplate surface has an evenly distributed heater element providing homogeneous heat distribution.

The non warping AISI 304 stainless steel sand bath can be easily removed for cleaning.

The Hotplate can be used independently without the sand bath. The hotplate casing is made of AISI 304 stainless steel.

### CONTROL PANEL

Mains power switch with luminous “On” indicator lamp.

Electronic heater power control between 0 to 100 %.

Heater operation indicator lamp.

The back of the unit is fitted for an electronic thermometer “Sensotherm” or digital electronic controller “Electemp-TFT” described under regulators and controllers (see page 293).



### MODELS

Part No.	Height / Width / Depth (usable) cm	Height / Width / Depth (exterior) cm	Power W	Weight Kg
6000708	5 18 38	18 22 41	2300	13
6000709	5 30 30	21 32 35	2800	15



## Circulation pump “Ciragit”

### FEATURES

The unit can be adjusted to fit most tanks or baths, can be submerged and held in position by an adjustable height support clamp.

Circulation pump allowing external circulation of liquids. ON / OFF illuminated power switch with light when on.

### SPECIFICATIONS

OUTPUT water:

20 °C, 12 l/ min.- 150 mbar

OUTPUT silicon oil:

50 °C 4 l/min.- 130 mbar

Dimensions:

9.5 x 12.5 x 23 cm

Power: 30 W

Weight: 2.5 Kg.

Part No. 3000402



NEW DESIGN