



Microscope slide dispenser

Code - BPG012

Nahita's microscope slide dispenser has been specifically developed for you to store and preserve your slides from any damage.

Thanks to this little dispenser you will easily withdraw one slide at a time just by rotating the knobs that are located at both sides of this item. Thus, you will obtain dust-free microscope slides without putting your fingerprints on them.

Depending on the thickness of the slides loaded the dispenser can store up to 80 slides. As it has a transpa-



Technical specifications

Code	BPG012
Capacity	80 standard microscope
Dimensions	120x100x165 mm
Made of	ABS plastic