



Pasteur pipettes

An indispensable element in daily laboratory work for carrying out routine and experimental techniques that require the measurement of imprecise volumes, vacuum aspiration and even sowing on bacteriological plates. Our catalogue offers non-sterile Pasteur pipettes in both glass and polyethylene.

Glass Pasteur pipettes

Code **DLB001, DLB002, DLB020, DLB021**

- Made of soda-lime glass, with a long, fine tip, and with a throttle on the suction tube.
- Non-sterile
- Non graduated
- They have a capacity of approx. 2 mL.
- Ref.: DLB020 and BLB021 include a cotton filter at the upper end.



Technical specifications

Code	Length (mm)	Units / box
DLB001	150	250
DLB002	230	250
DLB020	150	250
DLB021	230	250

Teats for Pasteur Pipettes

Code **MBL013, MBL014, MBL015**



Technical specifications

Code	Model	Material
MBL013	Dummy	Latex
MBL014	Normal	Red rubber
MBL015	Strong	Red rubber