


## Volumetric Pipettes 1 mark, Class AS

- Compliance with USP standard
- With Batch Certificate

Item Code	Capacity (ml)	Colour Code	Tolerance ( $\pm$ ml)	Pack of
123.461.00	0.5	2 BLACK	0.006	5
123.461.01	1	BLUE	0.006	5
123.461.02	2	ORANGE	0.006	5
123.461.02A	3	BLACK	0.010	5
123.461.02B	4	2 RED	0.010	5
123.461.03	5	WHITE	0.010	5
123.461.03A	6	2 ORANGE	0.010	5
123.461.03B	7	2 GREEN	0.010	5
123.461.03C	8	BLUE	0.020	5
123.461.03D	9	BLACK	0.020	5
123.461.04	10	RED	0.020	5
123.461.04E	15	GREEN	0.030	5
123.461.05	20	YELLOW	0.030	5
123.461.06	25	BLUE	0.030	5
123.461.06E	30	BLACK	0.030	5
123.461.06O	40	WHITE	0.050	5
123.461.07	50	RED	0.050	5
123.461.08	100	YELLOW	0.080	5



ENQUIRE



## Pipette Volumetric, Class-A, Serialized

- Made from ASTM E-438, TYPE-1, BORO 3.3 GLASS
- Complies with ASTM E-969 standard
- Printed in amber colour
- With individual work certificate

Item Code	Capacity (ml)	Tolerance ( $\pm$ ml)	Colour Coding Band	Pack of
123.454.00D	0.5	0.006	2 BLACK	6
123.454.01	1	0.006	BLUE	6
123.454.02	2	0.006	ORANGE	6
123.454.02A	3	0.01	BLACK	6
123.454.02B	4	0.01	2 RED	6
123.454.03	5	0.01	WHITE	6
123.454.03A	6	0.01	2 ORANGE	6
123.454.03B	7	0.01	2 GREEN	6
123.454.03C	8	0.02	BLUE	6
123.454.03D	9	0.02	BLACK	6
123.454.04	10	0.02	RED	6
123.454.04E	15	0.03	GREEN	6
123.454.05	20	0.03	YELLOW	6
123.454.06	25	0.03	BLUE	6
123.454.06E	30	0.03	BLACK	6
123.454.06O	40	0.05	WHITE	6
123.454.07	50	0.05	RED	6
123.454.08	100	0.08	YELLOW	2



ENQUIRE

