



Cell strainers

Code PZH001, PZH002, PZH003

- For preparation of stem cells and tissue-derived primary cells
- Sterile, DNase/RNase and pyrogen free
- Polypropylene shell and nylon screen
- Suitable for 50mL centrifuge tubes Corning, Thermo, etc.



Technical specifications

Code	PZH001	PZH002	PZH003
Specification	100 μ m, 150 mesh	70 μ m, 220 mesh	40 μ m, 330 mesh
Color	Green	Orange	Purple
Packaging	100 pcs/carton	100 pcs/carton	100 pcs/carton