




Several Instruments



DIAMOND POINT PENCIL

For glass engraving.
Part No. **1001130**



LABORATORY VIBRO-ENGRAVER SIGNOGRAPH-COM 25 DIAMOND

For engraving all types of glass, metal, plastic material, ceramic etc. Mains voltage 220-230 V. 50Hz. Power consumption: 25 W. Comes complete with a set of 5 engraving tips: Diamond, carbide, tempered steel and two for decorative work.
Part No. **1001129**



CORK BORER SHARPENER

With a wooden handle, blade of tempered steel.
Part No. **7000004**

CORK BORER

Electrobrite finish.
Universal handle.
Part No.

7000846 6 bores.
7000848 8 bores.
7008412 12 bores.

Electrobrite finish.
Individual handle.
Part No.

7000856 6 bores.
7000858 8 bores.
7008512 12 bores.



TUNGSTEN SCRIBE

For glass engraving, slides, ceramics, plastic etc. Retractable tip and pocket clip.
Part No. **1001674**



WATER JET VACUUM PUMPS

Part No.	Material	Pressure	Pump rate	Aspiration
4000633	PVC	2 Kg	8 l/min	2,6 l/min
7000293	Metal	2 Kg	5 l/min	1 l/min
7000292	Metal	2 Kg	5 l/min	0,5 l/min



MINI CORK BORER

1001332 6 bores: 4, 4,9, 5,5, 6,4, 8 and 9,5 mm Ø.



GLASS TUBE CUTTER

For tubes with diameters up to 30 mm.
Part No. **1000016**



CORK PRESS

Model Crocodile shape.
Part No. **6000070**



BLOW PIPE

Student model.
Total length: 250 mm.
Part No. **1000595**




PIPETTE FILLER, RUBBER BULB.

With three glass ball valves for filling, fine adjustment and empty. Accepts pipettes from 6 to 10 mm Ø.
Part No. **1001733**



MORTARS WITH PESTLE

Made of grey steel. Only suitable for dry products.
Part No. Ø Usable Height
1004713 130 mm 65 mm
1004715 150 mm 75 mm



RICHARDSON RUBBER BULB

With net covered decompression bulb.
For blowing out bench top automatic burettes.
Part No. **1001168**