

Extraction Apparatus



Condenser, Dimroth for Soxhlet Apparatus

- Made from ASTM E-438, Type-1, Class-A boro 3.3 Glass
- Condenser with glass screw thread connector

Item Code	Cone	Effective Length (mm)	Extractor (ml)	Connector	Pack of
192.202.01	45/40	300	100	Glass	1
192.202.01A	45/40	400	250	Glass	1
192.202.02	60/46	400	500	Glass	1
192.202.03	71/55	400	1000	Glass	1


5.1
BORO
GLASS

WHITE
PRINTING

Chemical
Resistant



ENQUIRE



Condenser, Dimroth for Soxhlet Apparatus

- Made from ASTM E-438, Type-1, Class-A boro 3.3 Glass
- Condenser with plastic screw thread connector

Item Code	Cone	Effective Length (mm)	Extractor (ml)	Connector	Pack of
192.181.01	45/40	300	100	Screw Thread	1
192.181.01A	45/40	400	250	Screw Thread	1
192.181.02	60/46	400	500	Screw Thread	1
192.181.03	71/55	400	1000	Screw Thread	1

3.3
BORO
GLASS

WHITE
PRINTING

Chemical
Resistant

SCREW
CAP



ENQUIRE