



Stopcock made of teflon



BORO
3.3



BORO
3.3

Features

- 1 | Endo Glassware dropping funnels are cylindrical and can be used as separating funnels or for adding approximate volumes to solutions.
- 2 | For this purpose, they present a lower standard ground cone to be attached to flasks or narrow neck bottles.
- 3 | The stopcock is made of teflon resistant to acids and bases.

Technical features

| Code | Graduation | Capacity (approx.) | Ground |
|--------|------------|--------------------|--------|
| PQL001 | No | 25 mL | 14/23 |
| PQL002 | Yes | 250 mL | 29/32 |