

► 200 series

1 | The educational stereomicroscopes are ideal to use in teaching for the mineral or little animal observation that sure will delight to the most curious pupils

2 | With head inclined 45° and equipped with incident and transmitted light.



01 | Binocular head



02 | Bilateral controls for focusing



03 | Stage with slide clips.

Characteristics

- 1 | They have binocular head with adjustable interpupillary distance and dioptic correction.
- 2 | Include a pair of eyepiece WF/10x
- 3 | Bilateral controls for focusing through rack and pinion mechanism.
- 4 | Black and white contrast stage and ground glass stage with slide clips.
- 5 | Objectives 2x and 4x that together with the eyepieces provide the necessary magnifications for the mineral and little insect visualization in great detail.

Code	HDB002
Model	218/2
Head	Inclined 45°
Eyepieces	WF10x/20mm (ø 30.5 mm)
Nosepiece	By pinion & rack
Interpupilar distance	55-75mm
Diopter correction	± 5 mm diopter
Objectives	2x/4x
Working distance	57 mm
Stage	B/N and ground glass
Ø Stage	95 mm
Illumination	Incident & trasmitted LED