



Square basin - 1 L - white

Ref. : G605621



Plus point

Basin fit for food contact.
Highly durable and chemical resistant.
Easy to clean.
Nesting with no suction effect.
Moulded rim makes it easy to grasp.

Product features

Customs code	39241000
EAN code	3100416056213
Net weight	0,09 kg
Length	160 mm
Width	160 mm
Height	68 mm
Capacity	1 L
Nesting rate	65,00 %
Product bottom internal height	120 mm
Product bottom internal width	120 mm
Product bottom internal height	65 mm
Material	LDPE
Colour	White
French made	Yes
Without Bisphenol A	Yes

Logistic data

Number of pieces per pallet	2160
Pallet weight	206,4 kg
Palette length	1200 mm
Pallet width	1000 mm
Pallet height	1760 mm
Number of units per layer	3
Number of layers per pallet	4
Number of pieces per sub-package	180
Weight in packaging	16,2 kg
Sub-package length	700 mm
Width under packaging	500 mm
Height under packaging	400 mm

All technical information on temperatures of use are given as an indication for the products concerned. The use of GILAC products remains the responsibility of the user.



Square basin - 3 L - white

Ref. : G605821



Plus point

Basin fit for food contact.
Highly durable and chemical resistant.
Easy to clean.
Nesting with no suction effect.
Moulded rim makes it easy to grasp.

Product features

Customs code	39241000
EAN code	3100416058217
Net weight	0,21 kg
Length	220 mm
Width	220 mm
Height	95 mm
Capacity	3 L
Nesting rate	74,00 %

Product bottom internal height	160 mm
Product bottom internal width	160 mm
Product bottom internal height	90 mm
Material	LDPE
Colour	White
French made	Yes
Without Bisphenol A	Yes

Logistic data

Number of pieces per pallet	1200
Pallet weight	264 kg
Palette length	1200 mm
Pallet width	1000 mm
Pallet height	1710 mm
Number of units per layer	20
Number of layers per pallet	6
Number of pieces per sub-package	10
Weight in packaging	2,1 kg
Sub-package length	220 mm
Width under packaging	220 mm
Height under packaging	330 mm

All technical information on temperatures of use are given as an indication for the products concerned. The use of GILAC products remains the responsibility of the user.



Square basin - 4.5 L - white

Ref. : G605921



Plus point

Basin fit for food contact.
Highly durable and chemical resistant.
Easy to clean.
Nesting with no suction effect.
Moulded rim makes it easy to grasp.

Product features

Customs code	39241000
EAN code	3100416059214
Net weight	0,28 kg
Length	260 mm
Width	260 mm
Height	110 mm
Capacity	4,5 L
Nesting rate	74,00 %

Product bottom internal height	195 mm
Product bottom internal width	195 mm
Product bottom internal height	105 mm
Material	LDPE
Colour	White
French made	Yes
Without Bisphenol A	Yes

Logistic data

Number of pieces per pallet	640
Pallet weight	191,2 kg
Palette length	1200 mm
Pallet width	1000 mm
Pallet height	1270 mm
Number of units per layer	16
Number of layers per pallet	2
Number of pieces per sub-package	20
Weight in packaging	5,6 kg
Sub-package length	260 mm
Width under packaging	260 mm
Height under packaging	610 mm

All technical information on temperatures of use are given as an indication for the products concerned. The use of GILAC products remains the responsibility of the user.



Square basin - 7 L - white

Ref. : G606021



Plus point

Basin fit for food contact.
Highly durable and chemical resistant.
Easy to clean.
Nesting with no suction effect.
Moulded rim makes it easy to grasp.

Product features

Customs code	39241000
EAN code	3100416060210
Net weight	0,49 kg
Length	300 mm
Width	300 mm
Height	130 mm
Capacity	7 L
Nesting rate	74,00 %
Product bottom internal height	225 mm
Product bottom internal width	225 mm
Product bottom internal height	125 mm
Material	LDPE
Colour	White
French made	Yes
Without Bisphenol A	Yes

Logistic data

Number of pieces per pallet	540
Pallet weight	276,6 kg
Palette length	1200 mm
Pallet width	1000 mm
Pallet height	1635 mm
Number of units per layer	12
Number of layers per pallet	3
Number of pieces per sub-package	15
Weight in packaging	7,35 kg
Sub-package length	300 mm
Width under packaging	300 mm
Height under packaging	575 mm

All technical information on temperatures of use are given as an indication for the products concerned. The use of GILAC products remains the responsibility of the user.



Square basin - 10 L - white

Ref. : G606121



Plus point

Basin fit for food contact.
Highly durable and chemical resistant.
Easy to clean.
Nesting with no suction effect.
Moulded rim makes it easy to grasp.

Product features

Customs code	39241000
EAN code	3100416061217
Net weight	0,62 kg
Length	350 mm
Width	350 mm
Height	135 mm
Capacity	10 L
Nesting rate	74,00 %
Product bottom internal height	270 mm
Product bottom internal width	270 mm
Product bottom internal height	130 mm
Material	LDPE
Colour	White
French made	Yes
Without Bisphenol A	Yes

Logistic data

Number of pieces per pallet	360
Pallet weight	232,2 kg
Palette length	1200 mm
Pallet width	1000 mm
Pallet height	1525 mm
Number of units per layer	9
Number of layers per pallet	4
Number of pieces per sub-package	10
Weight in packaging	6,2 kg
Sub-package length	350 mm
Width under packaging	350 mm
Height under packaging	430 mm

All technical information on temperatures of use are given as an indication for the products concerned. The use of GILAC products remains the responsibility of the user.



Square basin - 12 L - white

Ref. : G606221



Plus point

Basin fit for food contact.
Highly durable and chemical resistant.
Easy to clean.
Nesting with no suction effect.
Moulded rim makes it easy to grasp.

Product features

Customs code	39241000
EAN code	3100416062214
Net weight	0,71 kg
Length	390 mm
Width	390 mm
Height	140 mm
Capacity	12 L
Nesting rate	75,00 %
Product bottom internal height	290 mm
Product bottom internal width	290 mm
Product bottom internal height	135 mm
Material	LDPE
Colour	White
French made	Yes
Without Bisphenol A	Yes

Logistic data

Number of pieces per pallet	240
Pallet weight	179,4 kg
Palette length	1200 mm
Pallet width	800 mm
Pallet height	1570 mm
Number of units per layer	6
Number of layers per pallet	4
Number of pieces per sub-package	10
Weight in packaging	7,1 kg
Sub-package length	390 mm
Width under packaging	390 mm
Height under packaging	445 mm

All technical information on temperatures of use are given as an indication for the products concerned. The use of GILAC products remains the responsibility of the user.