

Bottle



Bottle, Gas Washing With Interchangeable Ground Joints

- Complies with DIN 12463
- Made from Heat & Chemical resistant Boro 3.3 Glass

Item Code	Capacity (ml)	Joint Size	Pack of
286.202.04	100	29/32	2
286.202.01B*	125	29/32	2
286.202.05	250	29/32	2
286.202.06	500	29/32	2

3.3
BORO
GLASS



*Non DIN



ENQUIRE



Gas Bottle Head

- Made from ISO/DIN 3585, Class A, TYPE 1, Boro 3.3 glass
- Complies DIN 12596 standard
- Made from Heat & Chemical resistant Boro 3.3 Glass

Item Code	Capacity (ml)	Joint Size	Pack of
285.202.04	100	29/32	1
285.202.01B*	125	29/32	1
285.202.05	250	29/32	1
285.202.06	500	29/32	1


3.3
BORO
GLASS



*Non DIN



ENQUIRE



Bottle, Gas Washing With Head

- Complies with DIN 12596 and 12463 standards

Item Code	Capacity (ml)	Joint Size	Pack of
285.406.01*	125	29/32	1
285.406.02	250	29/32	1
285.406.03	500	29/32	1

*Non standards



ENQUIRE