



Erlenmeyer flasks for cell culture

Code CDD007, CDD008, CDD009, CDD010

- Polycarbonate (PC) container, high-density polyethylene (HDPE) cap
- Electron beam sterilised
- Non-pyrogenic and DNase/RNase free
- With vent caps (0.2 μm filters)
- Individually wrapped in bags
- Flat bottom
- Can be used in a wide temperature range: -80°C - 121°C



Reference	Packaging
CDD007	1 pc/bag, 30 bags/carton
CDD008	1 pc/bag, 12 bags/carton
CDD009	1 pc/bag, 12 bags/carton
CDD010	1 pc/bag, 12 bags/carton

Technical specifications

Reference	CDD007	CDD008	CDD009	CDD010
Volume (mL)	125	250	500	1000