

Data Sheet

Product: Joseph paper 25 g/m²

Codes:

PSX001 Ream Joseph paper,15x15 cm, 500 sheets

PSX002 Ream Joseph paper, 35x50 cm, 500 sheets

- Composition: 100 % recycled cellulosic fibres
- Other components: -
- Method of manufacture: wet crepe, chlorine and acid free
- Color: white > 70 %
- pH: 7
- Recyclability rate: 100 %
- Certified FSC RECYCLE
- Grammage (NFQ 03-019): 25 g/m² ±5%

	Direction of motion	Cross direction
Breaking length (NFQ 03-004)	> 1000 m	> 900 m
Elongation (NFQ 03-004)	> 6.0 %	> 1.5 %

Signed:



AUXILAB, S.L.
N.I.F. ESB31072218
MATERIAL PARA LABORATORIO
Polig. Morea Norte, C/D. Huevo 6 - Tel. 34 948 310 513
31101 BERRIAIN-NAVARRA (ESPAÑA)

Technical Department