

Pharmacists' spatulas / Section lifters

SECTION LIFTER

Flexible blade, blade length 50 mm,
stainless steel 1.4310, with stainless steel handle



Total length in mm	150	150	150
Width in mm	10	15	20
	80.083.010	80.083.015	80.083.020

SECTION LIFTER

Flexible blade, blade length 50 mm,
stainless steel 1.4310, with wooden handle



Total length in mm	150	150	150
Width in mm	10	15	20
	80.063.010	80.063.015	80.063.020

PHARMACISTS' SPATULA

Extremely flexible blade,
chrome steel 1.4034, tempered



Blade length in mm	75	100	130	150
Blade width in mm	18	20	20	22
with wooden handle	80.861.075	80.861.100	80.861.130	80.861.150
with plastic handle, dishwasher safe	80.851.075	80.851.100	80.851.130	80.851.150

PHARMACISTS' SPATULA

Stainless steel blade, with wooden handle,
75 to 150 mm: Conical blade



Blade length in mm	75	100	130	150	200	250	300
Blade width in mm	18/22	18/22	18/22	18/22	30	30	30
Stainless steel 1.4310, flexible	80.061.075	80.061.100	80.061.130	80.061.150			
Chrome steel 1.4034, tempered	80.661.075	80.661.100	80.661.130	80.661.150	80.661.200	80.661.250	80.661.300

PHARMACISTS' SPATULA

Stainless steel blade, with plastic handle,
conical blade, dishwasher safe



Blade length in mm	75	100	130	150
Blade width in mm	18/22	18/22	18/22	18/22
Stainless steel 1.4310, flexible	80.051.075	80.051.100	80.051.130	80.051.150
Chrome steel 1.4034, tempered	80.651.075	80.651.100	80.651.130	80.651.150