



Infrared tip sterilizer "Sterilbio"

APPLICATIONS

Sterilization of culture tips, needles and small instruments. Incinerates organic material at 900 °C for 5/8 seconds.

FEATURES

Made from ceramic elements with a quality quartz tube interior.
Does not need purged oxygen and can be used in anaerobic environments.
Eliminates problems of aerosols and bacterial contamination.
Support at the side for 4 Kolle handles, needles, lancet or dissecting forceps etc.

MODEL

Part No.	Height/Width/Depth (exterior) cm			Ø/Length (tube interior) cm		Power W	Weight Kg
3000788	20	13	18	1.4	15	110	1.5



COMECTA Hot Air gun "Salki 1500"

WORKING TEMPERATURE 400 °C AND 590 °C.

APPLICATIONS

Suitable for heating liquids, solids and gasses:

- Dry
- Soft solder
- Fusions
- Evaporation
- Curing adhesives
- Plastic moulds
- Softening surfaces
- Fast defrost
- Film drying R.X.
- Chromatography applications drying paper and electrophoresis.
- Heating moulds and drying dental moulds etc.

FEATURES

Industrial model with vertical angle support.

Brushless motor.

Three positions:

- Position 0 stop.
- Position I temperature 400 °C.
- Position II temperature 590 °C.

MODEL

Part No.	Temperature °C	Power W	Weight Kg
5000750	until 590	1500	0.8



ACCESSORIES



Burner head 75 mm.
Part No. **5011148**



Shaft Ø 20 mm.
Part No. **5014415**



Tube deflector.
Surface 75 mm.
Part No. **5014449**

Diverse laboratory instruments

COMECTA Pocket Balance "Pocket Scale"

FEATURES

Capacity up to 200 g, measuring in g, oz and qty.

Digital LCD display.

Stabilization time ±3 seconds.

Tare by subtraction to maximum weight.

Piece counting function and automatic reset to zero.

Automatic calibration, external.

Batteries 3 x 1.5 V.

Lid to cover the weighing pan and operation menu.

MODEL

Part No.	Max. weight. capacity g	Precision ±g	Height / Width / Depth (exterior) mm			Width / Depth (plate) mm		Weight Kg
5830024	200	0.1	24	74	120	74	71	0.13

Comes with protective cover.

