

Funnel, Dropping, Cylindrical, Graduated With Glass Stopcock

- Made from DIN ISO 3585, BORO 3.3 GLASS
- With stem inside cone
- Complies with DIN 4800 Standards
- With PE Stopper

Item Code	Capacity (ml)	N/S	Bore (mm)	Stopcock Type	Pack of
165.210.02	50	14/23	2.5	Solid Glass	1
165.210.02A*	50	29/32	2.5	Solid Glass	1
165.210.03	100	14/23	2.5	Solid Glass	1
165.210.03A*	100	29/32	2.5	Solid Glass	1
165.210.05	250	29/32	4	Solid Glass	1
165.210.06	500	29/32	4	Solid Glass	1
165.210.07	1000	29/32	6	Solid Glass	1

3.3
BORO
GLASS

Chemical
Resistant

WHITE
PRINTING

GLASS

PE

WITH CONE
& SOCKET



ENQUIRE

* Non DIN



Funnel, Dropping, Cylindrical, Graduated With PTFE Stopcock

- Made from DIN ISO 3585, BORO 3.3 GLASS
- With stem inside cone
- Designed from DIN 4800 Standards
- With PE Stopper

Item Code	Capacity (ml)	N/S	Bore (mm)	Stopcock Type	Pack of
167.210.02	50	14/23	2.5	PTFE	1
167.210.02A*	50	29/32	2.5	PTFE	1
167.210.03	100	14/32	2.5	PTFE	1
167.210.03A*	100	29/32	2.5	PTFE	1
167.210.05	250	29/32	4	PTFE	1
167.210.06	500	29/32	4	PTFE	1
167.210.07	1000	29/32	6	PTFE	1

3.3
BORO
GLASS

Chemical
Resistant

WHITE
PRINTING

PTFE

PE

WITH CONE
& SOCKET



ENQUIRE

* Non DIN



Funnels, Pressure Equalizing, Pear Shape

- Made from DIN ISO 3585, BORO 3.3 GLASS
- With stem inside cone

Item Code	Capacity (ml)	N/S	Bore (mm)	Stopcock Type	Pack of
171.202.01	50	14/23	2.5	PTFE	1
171.202.02A	100	14/23	2.5	PTFE	1
171.202.04A	250	29/32	4	PTFE	1
171.202.05A	500	29/32	4	PTFE	1

3.3
BORO
GLASS

Chemical
Resistant

WHITE
PRINTING

PTFE

WITH CONE
& SOCKET



ENQUIRE

