




RUBBER STOPPERS WITH FOLDING SKIRT

Code - **BMJ012_ BMJ019**

Rubber stoppers with skirts are designed to ensure effective sealing in laboratory applications. They are made of flexible rubber, allowing them to adapt to different types of containers.

- Approximate working temperature range: 60–70 °C
- Manufactured in natural rubber
- Autoclavable on a one-off basis (once or twice at most).
Not suitable for repeated autoclave sterilization processes
- The material can be perforated without compromising its functionality



Code	Size* (mm)	D1 (mm)	D2 (mm)	H1 (mm)	H2 (mm)
BMJ012	10.2	10.2	13.3	9.7	10.2
BMJ013	13	13	18.7	13.5	13.4
BMJ014	15.85	15.85	20.2	15.2	15
BMJ015	24	24	28.3	13.5	21.4
BMJ016	13.6	13.6	18	10	15.6
BMJ017	18.5	18.5	24	14.8	18.5
BMJ018	24.4	24.4	28.8	15.4	21
BMJ019	28.6	28.6	34.65	19.75	16.5

*The size specification refers to the maximum inner diameter of the container to be sealed