

ENTONNOIRS CONIQUES AVEC POIGNEE

CONICAL FUNNELS WITH HANDLE

Destiné à transvider divers produits.

Extérieur bord roulé qui assure une meilleure solidité et évite les déformations.

Intérieur avec filtre non amovible, crochet à l'extérieur pour meilleure préhension et suspension.

Le tube d'écoulement est conique et présente une déformation dans la masse permettant une régularité dans l'écoulement, pour laisser l'air s'échapper. Il est étroit et peut être utilisé dans les bouteilles.

Intended to transfer various products.

Outside rolled edge that ensures a better solidity and avoids distortions.

Interior with a fixed filter, handle for better grip and suspension.

The flow tube is conical and has a deformation in the mass allowing a regularity in the flow, to let the air escape. It is narrow and can be used in bottles



Réf : 4215

Réf : 4220

Réf : 4250

Réf : 4215 Capacité 0 ,5 ml Ø 150 mm

Réf : 4220 Capacité 1L Ø 195 mm

Réf : 4250 Capacité 3L Ø 250 mm

Ref: 4215 Capacity 0.5 ml Ø 150 mm

Ref: 4220 Capacity 1L Ø 195 mm

Ref: 4250 Capacity 3L Ø 250 mm